

Padded Swivel Seat Cushion



HCPCS:



FIGURE A



FIGURE B

- Swivel seat allows individual to turn up to 360 degrees
- Makes getting on or off a seat or chair easy and comfortable (Figure A)
- Durable and stable
- Portable: weighs only 3.25 lbs
- Dimensions: 15 /" diameter
- Retail Packaged
- Weight Capacity 300 lbs

PRODUCTS

Description

RTLAF-300	Swivel Seat Cushion,Gray,RTL,1EA,8/CS
-----------	---------------------------------------

SPECIFICATIONS

Description

Warranty	Limited 1 Year
Weight Capacity	300 lbs
Weight Capacity (13046 13046-HD)	300 lbs.