

CHAD® 3-in-1 Oxygen Cylinder Shoulder Carry Bag

HCPCS:



- Hook and Loop Fastener brand touch storage pouch
- Opening for threading cannula on side of bag
- Vertical, horizontal, or backpack configuration
- Compatibility: M4, M6, M7, M9, while using any Drive or CHAD® conservers and regulators
- Size: 15.75" (H) x 4.33" (W) (at base) x 8.00" (W) (at widest point)

PRODUCTS

Description

OP-150-800

Carrying Bag, Vertical-Horizontal,1/ea