

CHAD® Regulators

HCPCS:



FIGURE A



FIGURE B

- 20-micron, sintered-bronze inlet filter
- FDA-recommended viton and brass seal washers
- All flows tested to accuracy within +/- 15% tolerances
- Meets ASTM-G175-03 promoted ignition tests
- Inlet pressure: 200 3,000 PSI
- Outlet pressure: 50 PSI
- 100% QC tested in the USA

PRODUCTS

Description

CH3115-L	RegCHAD,BrassCre,0-15L,CGA870,Barb,1/ea
CH4808-L-BLUE	RegCHAD,Blue,0-8L,CGA870,Barb,1/ea
CH4815-L-BLUE	RegCHAD,Blue,0-15L,CGA870,Barb,1/ea

SPECIFICATIONS

Description

Warranty	Limited 3 Year
----------	----------------

ACCESSORIES

Description

OM-900-FEC	E-Clip for Conserver Handles, 10/bag, CS
SEAL3000	Seal Washer, CHAD O2 Cnsrvrs&Regs,10/pkg, CS