

DeVilbiss iFill® Cylinders

HCPCS:



FIGURE A



FIGURE B

- These proprietary cylinders are for use with the DeVilbiss iFill® Personal Oxygen Station
- Available in various sizes with 870 Post Valve or with integrated regulators (Figure A-C)

PRODUCTS

Description

| | |
|-------------|--|
| 535D-C-870 | C CYL W/IFILL ADPT AND 870 POST VALVE |
| 535D-D-870 | D CYL W/IFILL ADPT AND 870 POST VALVE |
| 535D-D-CF | D CYLINDER W/CF REGULATOR |
| 535D-E-870 | E CYL W/IFILL ADAPTOR AND 870 POST VALVE |
| 535D-E-CF | E CYLINDER W/CF REGULATOR |
| 535D-M6-870 | M6 CYL W/IFILL ADPT AND 870 POST VALVE |

ACCESSORIES

Description

| | |
|---------|----------------------------|
| 9286-RD | SEAL REGULATOR, 10 PKG, CS |
|---------|----------------------------|