

## DeVilbiss iFill® Cylinders

HCPCS:



FIGURE A



FIGURE B

- These proprietary cylinders are for use with the DeVilbiss iFill® Personal Oxygen Station
- Available in various sizes with 870 Post Valve or with integrated regulators (Figure A-C)

### PRODUCTS

#### Description

535D-C-870	C CYL W/IFILL ADPT AND 870 POST VALVE
535D-D-870	D CYL W/IFILL ADPT AND 870 POST VALVE
535D-D-CF	D CYLINDER W/CF REGULATOR
535D-E-870	E CYL W/IFILL ADAPTOR AND 870 POST VALVE
535D-E-CF	E CYLINDER W/CF REGULATOR
535D-M6-870	M6 CYL W/IFILL ADPT AND 870 POST VALVE

### ACCESSORIES

#### Description

9286-RD	SEAL REGULATOR, 10 PKG, CS
---------	----------------------------