

Pediatric Fingertip Pulse Oximeter

HCPCS:



FIGURE A

- Four-way OLED with six display modes
- Prescription required
- Portable, non-invasive spot check of oxygen saturation and pulse rate
- Uses two AAA batteries
- Low battery indicator
- Comes in retail box
- Easy operation auto-off feature
- Lanyard, silicon case, and hang bag included
- Child-friendly design
- Adjustable brightness levels from 1 - 10
- Finger size range: 0.3" - 1.0"

PRODUCTS

Description

| | |
|-------|--------------------------------|
| 18707 | Pediatric Pulse Oximeter, 1/ea |
|-------|--------------------------------|

SPECIFICATIONS

Description

| | |
|----------|----------------|
| Warranty | Limited 2 Year |
|----------|----------------|