

## CHAD® Evolution Motion Auto-Adjusting Oxygen Conserver

HCPCS:



- Delivers a uniform oxygen pulse at increments: 1, 2, 3, 4, 5 and 6 Plus CF preset @ 2 LPM
- Senses patient activity and turns oxygen flow up or down as needed, within seconds
- Program a Rest and Active setting equivalent to 1-6 liters per minute
- Minimum 1-year battery life\* utilizing 2 AA batteries
- Delivers a 5:1 savings ratio from 14-40 breaths per minute
- Senses a breath/delivers oxygen in the first half of the inspiratory cycle
- Operating pressure range: 200-3000 PSI
- Weighs only 14.9 ounces with batteries
- Single push button selects the equivalent LPM dosage settings
- Auto shut-off feature conserves battery life
- Fits easily in virtually any carry bag

### PRODUCTS

### Description

OM-900M	Evolution Motion Conserver, 1/ea
---------	----------------------------------

### SPECIFICATIONS

### Description

Warranty	Limited 2 Year
----------	----------------

### ACCESSORIES

### Description

OM-900-FH	S-Handle for OM-900/800/500 Series, 1/e, CS
OM-900-BD	Battery Doors & Screws, OM-900/OM-900M, 5/pkg, CS
OM-900-FEC	E-Clip for Conserver Handles, 10/bag, CS
SEAL100	Seal Washer, Conserver/Regulators, 10/pkg, CS