

Extra Soft Cozy Cannula

HCPCS:



- Designed with curved prongs and an extremely soft cannula with a variety of tubing lengths to accommodate a wide number of patients' needs
- Compatible with supplemental oxygen sources delivering up to 15 LPM
- Latex free

| PRODUCTS | Description |
|------------|----------------------------------|
| SOFT 200 | Ad COZY cannula only CS/50 |
| SOFT 204 | COZYcannula w/4tub CS/50 |
| SOFT 207 | COZYcannula w/7tub CS/50 |
| SOFT 207 P | Pedi COZYcannula w/7tub CS/50 |
| SOFT 214 | Ad Cozy cannula w/14' tub, 50/cs |
| SOFT 225 | Ad COZYcannula w/25' tub 25/cs |
| SOFT 250 | Ad COZY cannula w/50 tub CS/20 |