

iGO2 Portable Oxygen Concentrator (POC)

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G

- Settings 1-5 Pulse Dose
- SmartDose Auto-Adjusting Converter Technology included
- Rugged design, including protective overmold
- Weighs less than 5 lbs
- Most sensitive converter trigger on the POC market
- 8.4" (H) x 3.5"(W) x 8.6"(L)
- Portable Oxygen Concentrator iGo family

PRODUCTS

Description

125D	iGO2 Portable oxygen concentrator NA
125D-XB	iGO2 Portable oxygen conc w/xtra bat NA

SPECIFICATIONS

Description

Battery Duration	3.5 hours @ setting 2
HCPCS	E1390 + E1392
Oxygen Purity	90% +4/-3 over entire operating range
Sound Level	37db

ACCESSORIES

Description

DV51D-606	US POWER CORD, EA
DV6X-619	DC car adaptor for IntelliPAP 2, EA
125D-613	iGO2 Battery pack, EA
125D-670	iGO2 Carrying case, EA
125CH-613	iGO2 Battery Charger (US), EA