

iGO2 Portable Oxygen Concentrator (POC)

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G

- Settings 1-5 Pulse Dose
- SmartDose Auto-Adjusting Conserver Technology included
- Rugged design, including protective overmold
- Weighs less than 5 lbs
- Most sensitive conserver trigger on the POC market
- 8.4" (H) x 3.5"(W) x 8.6"(L)
- Portable Oxygen Concentrator iGo family

PRODUCTS	Description
125D	iGO2 Portable oxygen concentrator NA
125D-XB	iGO2 Portable oxygen conc w/xtra bat NAc
SPECIFICATIONS	Description
Battery Duration	3.5 hours @ setting 2
HCPCS	E1390 + E1392
Oxygen Purity	90% +4/-3 over entire operating range
Sound Level	37db
Warranty	Drive DeVilbiss warrants this equipment to be free from defects in workmanship and materials for five (5) years from date of factory shipment to the original purchaser (typically the healthcare provider) with the exceptions of the Compressor (warranted for 3 years), Sieve Beds (warranted for 1 year), and Battery (warranted for 1 year). The carry case is excluded from this warranty. These warranties are limited to the Buyer of new equipment purchased directly from Drive DeVilbiss, or one of its Providers, Distributors, or Agents. Drive DeVilbiss obligation under this warranty is limited to product repair (parts and labor) at its factory or at an Authorized Service Center. Routine maintenance items, such as filters, are not covered under this warranty, nor does it cover normal wear and tear.

ACCESSORIES	Description
DV51D-606	US POWER CORD, EA
DV6X-619	DC car adaptor for IntelliPAP 2, EA
125D-613	iGO2 Battery pack, EA
125D-670	iGO2 Carrying case, EA
125CH-613	iGO2 Battery Charger (US), EA
DV68-620	iGO2 AC power supply, EA