

## Oxygen Cylinder Metal Stand

HCPCS:




FIGURE A

- Reversible clamp adjusts to hold cylinders from 7.25" to 9" in diameter
- Knock-down construction allows for a smaller package size, ideal for shipping and stocking
- Sturdy metal construction
- For cylinder sizes: M60, M90, MM, H & T

### PRODUCTS

|         | Description                             |
|---------|---|
| 18149KD | Cylinder Stand (KD) M60, M90,M120, 1/ea |

### SPECIFICATIONS

|           | Description         |
|-----------|---------------------|
| Dimension | Boxed   Assembled   |
| Height    | 5.1"   11.4"        |
| Length    | 19.7"   18.1"       |
| Warranty  | Limited Lifetime    |
| Weight    | 12.1 lbs   8.0 lbs. |
| Width     | 19.2"   17.0"       |