

Suction Therapy Accessories

HCPCS:

| PRODUCTS | Description |
|--------------|--|
| 4650D-609 | PUSH-IN FEET (4-PK) |
| 7305D-605 | SEP-SALES COMPRESSOR ASS'Y |
| 7305D-606 | SEP SALES CARRYING CASE |
| 7305D-607 | REPLACEMENT GAUGE |
| 7305D-608 | BACTERIA FILTER (NON-STERILE) 12EA/PKG |
| 7305D-609 | SEP-SALES HEAD CAP ASS'Y |
| 7305D-610 | CIRCUIT BOARD, 7305D SERIES |
| 7305D-612 | REPLACEMENT TUBING 6-PK |
| 7305D-614 | ElbowConct,Filter, 3.375" tubing, 1ea/ea |
| 7305d-615 | CUP SEAL (DFT) SEPARATE SALE PKG |
| 7305d-616 | BOT COVER ASSY (DFT) SEP SALE PACK |
| 7305D-617 | TOP COVER SEPARATE SALE PACKAGE |
| 7305d-618 | LED LABEL -DFT SEPARATE SALE PACKAGE |
| 7305D-619 | ELBOW (DFT) 6-PK SS |
| 7305D-623 | CYLINDER W/VALVES |
| 7305D-632 | Int Fltr Disp Canstr PK48 |
| 7305D-633 | Int Filter Disp Canister Comp Kit Assmd |
| 7305D-633-25 | Int Fltr Disp Cnstr Comp Kit Assmd Cs 25 |
| 7305D-635 | 800ML DISP CONTR FILTER Cart. Asml,12 P |
| 7305P-605 | LED LABEL 5 PK |
| 7305p-610 | CIRCUIT BOARD, 7305P SERIES |
| 7305P-614 | SEP SALES BATTERY ASSY,1/ea |
| 7305p-620 | BATTERY RETROFIT KIT |
| 7305P-631 | CONTINENTAL EUROPE CORD |