

## IntelliPAP® CPAP Pillow

HCPCS:



FIGURE A



FIGURE B

- Intelligent design offers space for mask and tubing to reduce mask drag and shift, as well as tubing interference while sleeping
- Provides better head and neck comfort while sleeping with a CPAP mask
- Comes with machine washable pillow case
- Multiple positioning options increases comfort for back, side and stomach sleepers
- Made in the USA of hypoallergenic materials

### PRODUCTS

#### Description

CPAP-PILLOW	CPAP Pillow 1/ea, 6/cv
CPAP-PILLOWCASE	Replacement Pillow Case 1/ea, 6/cv

### SPECIFICATIONS

#### Description

Dimensions	26" x 16" x 5"
------------	----------------