

LS9500 10" Lateral Rotation Mattress with on Demand Low Air Loss

HCPCS: E0277



- Combination therapy active system incorporates the functions of lateral rotation and low air loss to assist with pulmonary blood flow and mucosal drainage as well as providing microclimate control to keep the patient cool and dry
- Designed to prevent, treat and heal pressure ulcers in the home or long term care setting
- Control unit offers 10 comfort levels, 10,20, 30 & 60 minute adjustable turning cycles and maximum turn angle adjustments
- Turning modes include combinations of ¼, ½, ¾, and adjustable turning positions up to 40 degrees with left, right, both or none (static) directions
- Audio and Visual alarm alerts when air pressure is low
- 10" high side air bolsters provide fall prevention protection and facilitates safe patient transfers
- Pillow feature provides head stability and comfort
- Advanced Clinical and Care Giver Features: Auto-Firm, auto-recovery, microprocessing feedback, lock-out and more!
- Unique split air cells provide customizable turn time and rotation angles for patient comfort and satisfaction
- Solenoid valves prevent air loss and provide power outage protection
- Removable quilted Nylon cover is fluid resistant, low shear and vapor permeable helping protect the skin from friction and moisture
- CPR valve allows for rapid deflation

PRODUCTS	Description
LS9500	LateralRotationMattressw/LAL, 1/ea

SPECIFICATIONS	Description
Control Unit Airflow	45 LPM
Control Unit Alarms	Visual/Audible
Control Unit Type	Compressor
Mattress Cover Material	Fluid-Resistant Nylon
Warranty	1 Year
Weight Capacity	500 lbs

ACCESSORIES	Description
LS9525	Full Mattress Replacement, LS9500, 1/ea, EA
LS9520	Pump Only, LS9500, 1/ea, EA
LS9540	Replacement Cover, LS9500, 1/ea, EA
LS9550-CELL	Individual Bladder Only, LS9500, 1/ea, EA