

Multi-Ply ShearCare 1500 Pressure Redistribution Foam Mattress

HCPCS: E0184



FIGURE A

- Bariatric, latex free, therapeutic foam mattress provides pressure redistribution and shear/friction reduction
- Full width top layer is high-resilient, high-density die-cut foam for pressure redistribution and maximum support
- High-density foam bottom layer provides support and stability for the bariatric patient
- Patented LuraQuilt foam heel section offloads pressure from vulnerable heels
- Removable premium, multi-stretch cover is fluid resistant, low shear and vapor permeable helping protect the skin from friction and moisture
- Meets flammability standard CFR 1632 and CFR 1633
- Shear Care 1500 Pressure Redistribution Foam Mattress

| PRODUCTS | Description |
|--------------|---------------------------------------|
| 1500SC-42-FB | ShearCare 1500 Mattress 42x80, 1/ea |
| 1500SC-48-FB | ShearCare 1500 Mattress 48x 80", 1/ea |

| SPECIFICATIONS | Description |
|------------------------------------|---|
| Cover Material | Fluid Resistant Nylon Top Waterproof Vinyl Bottom |
| Dedicated Heel Section | Yes |
| Mattress Layers of Foam | 2 |
| Side Rails | No |
| Warranty | 5 Year Limited Foam 3 Year Cover |
| Weight Capacity (13046 13046-HD) | 42" 650 lbs, 48" 750 lbs |

| ACCESSORIES | Description |
|-------------|---|
| CPU64280 | Cover,MultiPly/ShearCare1500Ser,6x42x80, EA |
| CPU64880 | Cover,MultiPly/ShearCare1500Ser,6x48x80, EA |
| CPU65480 | Cover,MultiPly/ShearCare1500Ser,6x54x80, EA |