

Walker Accessory Bag

HCPCS:



- Adjustable bag loop fasteners attach to wide range of mobility aids
- Spacious interior pocket with zipper closure and key ring
- Three front pockets hold essentials such as keys, snacks or cell phone
- Two back pockets allow items to be accessed from users side
- Lined interior and pockets
- Treated with water resistant technology
- Spot clean with wet cloth
- Take your items with you wherever you go using the functional and versatile bag that easily secures to most walkers. The front and rear pockets are perfect for keys or a cell phone, and the large lined interior pocket with O-ring zipper closure holds all your essential items. The bag is constructed of lightweight nylon material that will remain sturdy even with extensive use.

PRODUCTS	Description
RTL10256NV	Univ MobilityBag, Med Flat, Navy, 1EA
RTL10256PK	Univ MobilityBag, Med Flat, Pink, 1EA
RTL10256TL	Univ MobilityBag, Med Flat, Teal, 1EA

SPECIFICATIONS	Description
Dimensions	9" H x 16" W x 2.5" D
Retail Box Dimensions	6.1" L x 1.4" W x 11" H
Weight	.5 lbs
Weight Capacity	6 lbs.